

## ANSI/AAMI

Compliant




Ready to Keep pace with life in any surgical environment.

www.egicd.com

### CONTROL

The EGIC steam sterilizer's PLC Siemens continuously monitor and control all system operations and functions.


All control configuring is performed through the easy - to - use Siemens touch screen display 7 inch, from which all sterilizer functions are controlled The date includes: Date, daily cycle number, starting time of each cycle, temperature selected Equipped with thermal Printer which provides an easy -to- read recorded data of all printed cycle data.





### **Built-in Steam Generator**

316L stainless kw electric integral 54-18 The EGSTER Steam sterilizer is equipped with steel steam generator includes a manual flush and drain package





### **Door/Doors**

Made from a single formed piece of 316L stainless steel Sliding power The doors are equipped with a mechanical locking mechanism The doors are sealed by steam activated door seal gasket constructed of special type of long-life rubber compound that needs no lubrication



# E G I C ANSI/AAMI

Compliant

Increased shelf clearance allows processing Larger orto and bariatric sets and large retractors Requires no additional space

Model	Size (capacity lit)	Туре	
EGSTER 40	160	prevac	
EGSTER 50	280	prevac	
EGSTER 60	400	prevac	
EGSTER 70	670	prevac	

### Customizable cycles/cycle names variety of per-programmed cycles

Available Cycles:	Sterilize Temp.	Sterilize Time.	Dry Time.
Immediate-use Prevac	°270F (°132C)	4Min	3Min
Prevac	°270F (°132C)	3Min	30Min
Prevac	°270F (°135C)	4Min	5Min
Prevac	°270F (°132C)	4Min	20Min
Prevac	°270F (°132C)	30Min	15Min
Gravity	°250F (°121C )	15Min	30Min
Gravity	°270F (°132C)	30Min	30Min
Gravity	°250F (°121C)	25Min	15Min
Liquid	°270F (°132C)	45Min	N/A
Leak <mark>Test</mark>	N/A	N/A	N/A
DART Test (Bowie Dick)	°270F (°132C)	3.5Min	1Min
DART Warm-up	°270F (°132C)	3Min	1Min





## Compliant. Larger capacity. Ready for immediate use

### Safe, intuitive operation

- 1-Water saving temperature control
- 2 Automatic Door/Doors
- 3-Colour touch screen with security access
- 4 Simple programming; customizable cycle names
- 5 Display messages in complete phrases
- 6 Thermal printer
- 7 Control at eye level for easy access
- 8 Environmental friendly vacuum pump

#### Adaptable Configurations

- Single door
- Double door
- Recesses
- Cabinet
- House steam

#### Reliably engineered

- Non-lubricated, steam activated door gasket
- Easy access modular piping
- Fully jacketed pressure vessel
- 316L stainless steel construction
- 20 year chamber warranty

#### Accessories

- Rack/Shelves
- Loading Car/ Transfer Carriage
- Water Softener